SCOOTER TRAVELER

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IMPORTANT

Read the instructions carefully before using this product and keep them if necessary.

User: If you are unable to read or understand any warnings, precautions or instructions, please contact your health care professional or distributor before using this product or risk and damage may result.

Consult your doctor or health care professional before using a walking aid. He/she will advise and inform you on the most appropriate walking technique for you.

Distributor: this leaflet must be given to the users of this product.

PRODUCT DESTINATION

The Traveler is a medical device. It is a 4-wheel scooter that can be dismantled without tools into 3 parts (including a quick-release battery) for easy access in a car boot. Comfort seat with fold-up, height-adjustable armrests. The scooter does not have handles but has oval directional handlebars.

PRODUCT DESCRIPTION

SAFETY PRECAUTIONS

Failure to follow the safety precautions may result in the risk of serious injury.

- Do not exceed the maximum weight limit of 120 kg.
- If any parts are damaged, broken or worn, contact your dealer.
- Please read this manual carefully before using the scooter.
- For your own safety and that of others, please observe local traffic regulations when using this product.
- Do not attempt to modify the scooter, as this may damage it. Any modification to the product will invalidate the warranty in the event of a problem.
- Please check the following safety points before using the scooter:
 - o Check the tyres for damage or excessive wear.
 - o Check all electrical connections and how they are connected.
 - o Check the brakes to ensure that they work properly.
 - o Check the battery level to avoid running out of power during use.
- Before using the scooter, the dealer must carry out a test of the user's driving ability.
- When using the scooter, do not exceed the maximum load. The seat can only support one person at a time. Failure to comply with this instruction may cause injury to people and to the vehicle.
- When using our product, do not attempt to drive over large obstacles. Keep upright when crossing obstacles.
- When using the scooter on uneven or wet ground, be sure to slow down. The scooter cannot float and is sensitive to water. If the water level rises as high as the battery, this could cause a short circuit and damage the electrical parts of your scooter.
- Please slow down when turning or driving on slopes. Be careful not to lean too far forwards or backwards when using the scooter, so as to maintain a constant driving speed.

INSTALLATION AND DISASSEMBLY

1. Components

2. Fitting the handlebars

- While holding the handlebars, unlock by turning the knob anti-clockwise and remount the handlebars, then lock the knob by turning it in the opposite direction.
- Firmly lock the knob

Front and rear assembly :

Aligning the chassis: Move the front section closer to the rear section until the hook on the front section comes into contact with the spar on the rear section and the lock engages automatically.

Fitting the seat tube :

Place the tube in the housing and lock it in place using the pin.

Fitting the battery pack :

Present the battery pack as shown below. The battery pack and power supply will connect automatically.

Fitting the seat :


- a. Place the seat on the seat tube and rotate it to fix its position.
- b. To fit the armrests, insert them into the tube under the seat and lock them in place using the dowel pins, then secure them using the tightening knobs.

3. Disassembly

Press the locking knob at the bottom of the handlebar and turn it to the right. See photo opposite: this will lock the handlebars.

- To remove the scooter seat, simply lift it up.
- To remove the batteries from the scooter, simply lift them up by the handle.


ADJUSTMENT AND USE

FUNCTIONAL ELEMENTS

A. Dashboard

The dashboard contains all the controls required to drive the scooter: the key switch, the travel buttons, the horn, the speed button and the battery indicator. See below.

WARNING! The dashboard must not be exposed to wet conditions. If the scooter is stored in a damp environment or in the rain, do not use it until the dashboard is completely dry.

1. Horn

Press this button to start the horn. To be used in the event of danger.

2. Key switch

Place the key in the key switch and turn clockwise to start the scooter. Turn the key anticlockwise to switch the scooter off.

CAUTION: If you turn the key to the "off" position while the scooter is running, the electronic brake will engage and the scooter will stop automatically.

3. Battery level meter

After starting the scooter, the battery status meter displays the battery power.

4. Setting the speed

You can use the speed control dial to select the scooter's minimum or maximum speed.

5. Travel buttons


You can use these buttons to control the forward and reverse direction of the scooter. Place your left and right hands on the left and right handle respectively.

Push the right button with the right thumb and the scooter will move forwards. Push the left button with the left thumb and the scooter will move backwards. When you release both hands from these buttons, the scooter will stop moving forward.

6. Light

You can use this button to control the scooter's lighting system. Push the button forwards and the light will come on. Push the button backwards and the light will go out.

B. Rear mechanical components

1. The Freewheel position lever

You can use this lever to switch the scooter to freewheel mode. Push this lever forward to stop the braking system. The scooter will then enter freewheel mode. You can now move the scooter with the engine stopped. Pull the lever back to reactivate the driving and braking systems. The scooter will no longer be in freewheel mode.

CAUTION Before putting the scooter into freewheel mode, be sure to remove the key from the switch. Do not sit on the scooter when using this lever. Do not place the scooter on slopes when in freewheel mode.

After using the freewheel mode, remember to return the scooter to the initial position (lever pulled back) to lock the brake.

2. Electromotor

The electromotor is a device that converts battery power into mechanical energy to drive the scooter's wheels.


3. Anti-rotation system

This system prevents the rear wheels from turning.


CAUTION! UNDER NO CIRCUMSTANCES SHOULD THIS SYSTEM BE REMOVED/DISASSEMBLED.

C. Battery pack

There is a handle at the top of the battery pack, making it easy to remove the battery from the scooter. The block contains two rechargeable batteries, the charger power input, the overload circuit breaker (reset button) and the battery fuse. See photo:

To recharge the battery, connect the charger to the scooter (using this power input) and plug it into a mains socket for 8 to 14 hours:

When the scooter's battery becomes weak, the main circuit breaker may trip to protect the motor and electronics from damage. If the circuit breaker trips, the entire electrical system will shut down and the reset button will appear where shown in the photo.


After tripping, allow the components to cool for 1-2 minutes, then press the reset button to restart the main circuit breaker. If the circuit-breaker trips frequently, please recharge the battery frequently, and contact the after-sales service if there is a battery problem.

Note: There is also a fuse in the dashboard to protect its electronic components. If one or more parts of the dashboard are not responding, check the fuse.

D. Transporting the scooter

When transporting the scooter in a vehicle, ensure that it is properly secured/attached. The photo below shows in black dots the parts that must be secured/attached/hooked up so that the scooter does not move when the transport vehicle is in motion.

CAUTION: Do not transport the scooter until it is securely fastened in the vehicle, as this could damage it.


COMFORT SETTINGS

1. Adjusting the height and tilt of the handlebars

To adjust the handlebars on your scooter, first remove the seat (by pulling it upwards) and the seat cover. Then adjust as shown in the photo:

Caution! Make sure you remove the key from the switch before adjusting the tilt. If your scooter is running, do not make any adjustments.


2. Seat orientation adjustment

Pull the rotating bar forward to release the seat.

- a. Turn the seat to the desired position.
- b. Return the seat's rotating bar to its initial position.

- ① Armrest adjustment knob
- ② Knobs
- ③ Rotating seat bar


3. Armrest adjustment

The armrests can be adjusted inwards or outwards.

- a. Lower the armrests to access the adjustment knob (point (1) in the photo above) and remove the knobs (point (2) in the photo above).
- b. Position the armrests inwards or outwards as required
- c. Replace the knobs (point (2) in the photo above) and raise the armrests (point (1) in the photo above).

4. Seat height adjustment

- ① Seat rod
- ② Lowest level of seat height
- ③ Position pin

- a. Remove the battery from the scooter, then the locking pin (point (3)).
- b. Raise or lower the seat to the required height.
- c. Replace the pin and battery.

USING THE BATTERY


1. Description

The traveler scooter is equipped with two long-life 12V 12 Ah batteries. To ensure safe and efficient charging, only use the charging system supplied with the scooter. Please do not charge the two batteries separately. The scooter should be charged for between 8 and 14 hours before it is used for the first time. To ensure stable and constant use: the scooter batteries must be fully charged.

Caution! Do not use an extension lead to connect the battery charger, use the standard power socket directly. Unattended children must not play near the scooter when the battery is being charged.

Do not charge the battery while the scooter is running. Keep the scooter away from sources of heat. Check that the charger, wiring and battery connectors are not damaged before using the scooter. If any part is damaged, contact the after-sales service. The temperature of the battery charger shell may rise during charging. To avoid burns, do not try to open it. Please do not expose the battery to adverse or extreme weather conditions. If it is inadvertently exposed to adverse conditions, set it aside before use.

The charge level of your batteries is indicated on the scooter's dashboard (from left to right): a green light indicates a full charge, a yellow light indicates a low charge and a red light means that the batteries are empty.

2. Frequently asked questions

a. How does my charger work?

When the scooter's charge is low, the charger works harder by transmitting more electric current to increase the amount of charge. When the battery is about to be fully charged, the charger will gradually reduce the amount of current it transmits. The charger's current is close to 0 amps when the battery is charged. The charger maintains the charge in the battery rather than overcharging the battery when connected.

b. Can another charger be used?

We supply a charger that exactly matches the type, capacity and properties of the battery pack. We therefore recommend that you only use the original charger supplied with the scooter. We would also remind you that it is not advisable to charge the two batteries separately.

c. What should I do if the scooter does not charge?

Check that the battery fuse is working correctly. Check that the plug is correctly connected to the scooter's power supply.

d. How long should the batteries charge?

There are two main factors to consider when charging the batteries.

- Daily use throughout the day
- Occasional use

Depending on how they are used, you will need more or less time to recharge the batteries. The battery charger we supply is designed not to overcharge. To ensure that the scooter works properly, please follow the instructions below:

- If you use the scooter every day, charge the batteries after each use. The batteries should be charged for 8 to 14 hours.
- If you use the scooter occasionally (once a week or once a month, for example), the batteries should be charged at least once a week for 24 hours.

Please note! Take care to keep the batteries fully charged and avoid overcharging. Refer to the instructions supplied with the charger.

If the scooter is used less than once a month, we recommend that you charge it for at least 48 hours to improve performance and ensure that the batteries last as long as possible.

e. How can battery life be extended?

Charging batteries to their full capacity maintains reliable capacity and extends their life. Keep the batteries fully charged as often as possible. Keep the scooter in a dry environment at room temperature. To prolong the life of the batteries, be sure to remove the key from the switch when the scooter is not in use.

f. How can I increase the number of kilometres per battery charge?

It is difficult to achieve ideal driving conditions: a smooth, hard, straight, flat road instead of a slope, pavements with cracks etc. Each condition will affect the maximum mileage per battery charge. Each condition will affect the maximum mileage per battery charge. Fully charge the battery before each use. Plan your route in advance if possible, to avoid obstacles in your path.

- Do not overcharge the scooter, only take the essentials.
- Try to maintain a constant speed.
- Avoid sudden stops and abrupt driving.
- Ensure that all cables are correctly connected.

g. What types of batteries are used, and what are their capacities?

We use two-cycle batteries. We recommend sealed, maintenance-free batteries such as AGM or colloidal batteries. Do not use batteries with removable covers.

Caution! Batteries contain corrosive chemicals. Please use AGM or colloidal batteries to reduce the risk of leakage or explosion.

h. Why does a new battery give a poor performance?

The chemical technology of batteries is different from that of car batteries or nickel cadmium cells. Deep cycle" batteries are specifically designed to provide fast power. (8 to 14 hours). New batteries are full when delivered, so they may experience temperature changes during transport. This may take a few days to stabilise.

i. Will the following steps ensure maximum efficiency and extended battery life?

1. Please fully charge the batteries before using them for the first time.
2. Use the scooter in a familiar and safe place. Drive slowly and do not go too far the first time.
3. Repeat steps 1 and 2 for 4 or 5 cycles for maximum efficiency and to prolong battery life.

j. Use on public transport?

Before using the scooter on public transport, please check that it is authorised.

Warning! Batteries contain lead and chemicals. Wear goggles and gloves if you have any problems with the batteries. Please wash your hands after handling.

Please note! The type, chemical properties, and ampere-hour (Ah) capacity of the two batteries must be exactly the same. Refer to the product specifications.

Keep tools and other metal objects away from the batteries. This can cause electric shock.

To replace the batteries, please use a spanner.

Caution! Do not attempt to lift weights beyond your physical reach. Please ask for assistance when setting up your scooter. We also advise you to ask your dealer for help with handling.

k. How do I change the batteries?

1. Remove the battery pack from the scooter and carefully place it on its side.
2. Remove the eight screws from the battery pack.
3. Place the positive battery up and remove the cover.
4. Disconnect the black (-) and red (+) electrodes from the batteries.
5. Remove the old batteries.
6. Place the new batteries in the battery pack.
7. Connect the red wire to the positive terminal of each battery (+).
8. Connect the black wire to the negative terminal of each battery (-).
9. Replace the cover and ensure that all the cables are securely inside the battery without being trapped.
10. Carefully turn the battery pack over.
11. Reinstall and tighten all the screws.

HOW IT WORKS

1. Before getting on the scooter

- Check that the batteries are correctly charged. *See Using the battery*
- Check that the "Freewheel" mode is not activated. *Cf. Rear mechanical components*
- Make sure that the handlebars are in the unlocked position: so that you can steer the scooter.

2. Getting on the scooter

- Make sure that the key is not in the switch. This will prevent your scooter from moving accidentally.
- Make sure the seat is securely fastened before you sit down.
- Make sure your feet are touching the base.

3. Adjustment and pre-drive check

- Is the seat at the right height?
- Is the seat properly secured?
- Is the key correctly inserted in the ignition? Is it turned to the "on" position?
- Does the horn work properly?
- Are there pedestrians, animals or obstacles in the way?
- Does your route include slopes or unstable terrain?

4. Driving

Be careful! When using the scooter, the following behaviours may affect your driving.

- Keep pets on a lead
- Carrying passengers or animals on the scooter
- Towing another vehicle

Warning! When using the scooter, always keep both hands on the handlebars and both feet on the base. This will enable you to drive your scooter at its best, and to stop in case of trouble.

- Set the ideal speed using the speedometer and lever - Parking brake
- To turn the scooter to the left, pull the handle to the left.
- To turn the scooter to the right, pull the handle to the right. If you need to stop, release the accelerator slowly. When the scooter has stopped, the brakes will automatically engage.

5. Getting off the scooter

- Make sure that the scooter has come to a complete stop before dismounting.
- Remove the key from the ignition.

6. Automatic stop function

In order to conserve the scooter's battery as long as possible, we have equipped our scooters with an automatic stop function. If you leave the key in the ignition in the "on" position but the scooter is not running, it will automatically switch off after 20 minutes.

If the automatic stop is activated, you can cancel it by following these steps:

- Remove the key from the ignition
- Reinsert the key to restart the scooter

TROUBLESHOOTING

Most scooters require regular checks. In general, most problems can be solved easily. Many problems are caused by a weak battery or damage to the batteries.

1. The alarm signal

The scooter's alarm signal helps you to quickly detect and rectify a fault. The alarm will sound in the following circumstances. The type of alarm varies from model to model

- If the alarm signal makes a series of fast and slow noises, this corresponds to a false signal. The alarm sounds only once. To stop the alarm signal, remove the key from the ignition, reinsert it and turn it to the "on" button.
- The alarm signal can identify a problem and sound it again after a short pause. The alarm will continue to sound until you switch the scooter off and on again.

Note: The scooter will not start until the fault has been rectified.

Below is a table of common errors you may encounter:

Alarm Signal	Problem	Solution
1. ■	Batteries very low, scooter will not start	Recharge as soon as possible
2. ■ ■	Low Batteries	Refill
3. ■ ■ ■	Voltage too high the scooter cannot be operated / Scooter is still connected to a socket outlet	Unplug the charger and restart the scooter. Be sure to reduce your speed when going downhill.
4. ■ ■ ■ ■	Protection activated	Switch the scooter off and restart it after a few minutes.
5. ■ ■ ■ ■ ■	The scooter is in freewheel mode	Switch off the scooter, remove the wheel mode free and restart the scooter
6. ■ ■ ■ ■ ■ ■	The accelerator lever is not in the middle when starting the scooter	Check that the battery is fully charged. If so, put the lever in the middle and restart the scooter.
7. ■ ■ ■ ■ ■ ■ ■ ■	Speed error	Contact your authorised distributor
8. ■ ■ ■ ■ ■ ■ ■ ■ ■ ■	Speed error	Contact your authorised distributor
9. ■ ■ ■ ■ ■ ■ ■ ■ ■ ■	Speed error	Contact your authorised distributor

2. What should I do if the scooter does not move when I press the accelerator lever?

If freewheel mode is activated to move the scooter manually: simply release the brake. To move ON the scooter, make sure freewheel mode is deactivated. Then try starting and accelerating again.

3. What should I do if all the scooter's systems stop responding?

- Make sure the key is in the ignition and in the "on" position.
- Check that the batteries are charged.
- Press the reset button on the main circuit breaker.
- Check that the batteries are correctly installed and that all the connections are secure.
- Check the dashboard fuse under the batteries.


4. Checking the fuse

Correct fuse / Blown fuse

- Remove the fuse from its slot
- Check whether the fuse is blown (see photo opposite)

- If the fuse is blown, fit a new fuse

- Caution! When replacing the fuse, make sure you replace it with a fuse of the same rating, otherwise it could damage the scooter.

5. What should I do if the scooter's batteries won't charge?

If the LED on the charger does not light up, check the battery fuse.

Replace the fuse if necessary.

- Make sure that both ends of the charger are correctly connected.
- If the fuse blows regularly, contact your dealer.

6. What happens if the battery power indicator does not respond when I start the engine or press the gear lever?

Charge the batteries fully

- Carry out a charge test by testing your batteries

If the scooter is not working properly and you cannot repair it, contact your supplier.

CARE AND MAINTENANCE

Failure to care for and monitor your product may result in the risk of serious injury.

Your scooter requires very little maintenance. The following parts require regular inspection and maintenance.

1. Tyres

Check tyres regularly for damage or wear.

Caution! Scooter wheels must be replaced by a qualified technician.

Caution! Make sure that the key is turned off and that the scooter is not in freewheel mode before carrying out the following steps.


2. Changing the wheel

The scooter is fitted with solid wheels. If a tyre is damaged, replace it with a new one.

3. Quick tyre repair

- Check that the scooter key is not in the ignition.
- Lift the side to be replaced. Put a piece of wood under the frame to raise the scooter.
- Remove the bolt and washer from the axle.
- Fit the new wheel
- Make sure that the key is in the correct position to avoid injury or damage to the scooter.

- Reinstall the bolt and washer. Tighten them securely

- ① Boulon - Screw
- ② Rondelle - Disc
- ③ Roue – Wheel
- ④ Clavette - Key
- ⑤ Rainure - Groove

4. Cleaning and disinfection

Use a damp cloth with a non-abrasive, non-corrosive product to clean plastic and metal parts.

Do not use cleaning products that could scratch the scooter.

Warning! Please observe the precautions for use of your cleaning product to avoid any risk to your skin or the scooter's paintwork.

5. Connections

Make sure that the connections are always well connected and that there is no corrosion.

6. ABS plastic

The scooter's bumper is made of very resistant plastic. You can polish it with car wax to keep it shiny.

7. Ball bearings and electric motor

These parts are pre-lubricated and sealed, so they require no maintenance.

8. Console, charger and rear electronics

These parts must be protected from moisture. If these parts are accidentally exposed to wet objects, do not use the scooter until they are completely dry.

9. Replacing the nylon lock nut

The old nylon locking nut must be replaced with a new one during maintenance, installation and de-installation. Nylon locknuts are not reusable as they are susceptible to damage, which reduces safety.

- Checklist:

To keep the product in good condition, check before each use, and check regularly as follows:

Check the product	Anytime	Weekly	Monthly
Check freewheel operation	X		
Checking your scooter's speed	X		
Scooter cleanliness			X
Check the condition of the battery			X
Check the condition of the tyre cable sheaths			X
Check that the seat is in the correct position for use	X		
Check the brakes	X		
Check the condition of the tyres		X	

CHANGE OF OWNERSHIP

The product can be reused. It must first be cleaned and disinfected in accordance with hygiene instructions. When transferring the product, remember to hand over all the necessary technical documents to the new user. Beforehand, the product must be inspected by an authorised specialist.

RECYCLING

When the product has become unusable and you have to throw it away: please contact your specialist dealer. If you wish to carry out recycling for yourself, ask your local waste management company for disposal guidelines.

STORAGE

Failure to observe the storage conditions may lead to deterioration of the product and thus to the risk of serious injury.

Do not store your product for prolonged periods near a source of heat or in direct sunlight (e.g. behind a window or near a radiator) or near a source of cold.

Keep away from open flames and sources of sparks.

Observe the following conditions when storing your rollator:

- In a dry, temperate place
- Protected by packaging, from dust and corrosion (e.g. abrasive elements, sand, seawater, salty air).
- Store all removed parts together, in the same place (or mark them if necessary) to avoid mixing them with other products when reassembling.
- All components must be stored without any load (do not place heavy objects on product parts, do not jam anything...).

If you are not going to be using the scooter for a while, we recommend the following:

- Fully charge the battery
- Remove the battery pack
- Store the scooter in a dry place away from extreme temperatures.
- If the battery is stored with a low charge, if you do not charge the battery often, or if you store the scooter in extreme temperatures, the battery may be seriously damaged.

Warning! Do not charge the battery if it has been stored below 0°C. This could cause serious damage to the battery.

When storing the scooter for a long period, it is advisable to protect the floor with a blanket or cardboard. This will reduce tyre wear.

TECHNICAL SPECIFICATIONS

Overall dimensions	Length 101 cm x Width 48 cm x Height 90 cm
180W motor 2 x 12V 12 Ah batteries Range 16km. Maximum speed 8km/h	
Weight	48 kg
Maximum load	120 kg

PERIOD OF USE

The life expectancy of this product is 5 years under normal use, safety and maintenance conditions. Beyond this period, the product can be used for as long as it is still in good condition.

WARRANTY

This product is guaranteed for two years from the date of purchase by the user. IDENTITIES will provide replacement parts for any part that has broken due to a manufacturing, design or material defect. Contact the distributor from whom you purchased the product in case you need to make a claim. Failures resulting from accidental damage, misuse, alteration, or wear and tear are not covered by this warranty. Contact the dealer from whom you purchased the product to apply the warranty.

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The Traveler is a medical device. It is a 4-wheel scooter that can be dismantled without tools into 3 parts (including a quick-release battery) for easy access in a car boot. Comfort seat with fold-up, height-adjustable armrests. The scooter does not have handles, but has oval directional handlebars.

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