**SCOOTER TRAVELER PLUS**

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## INSTRUCTIONS MANUAL EN

## IMPORTANT

Read the instructions carefully before using this product and keep them if necessary.  
User: If you are unable to read or understand any warnings, precautions or instructions, please contact your health care professional or distributor before using this product or risk and damage may result.  
Consult your doctor or health care professional before using a walking aid. He/she will advise and inform you on the most appropriate walking technique for you.  
Distributor: this leaflet must be given to the users of this product.

## PRODUCT DESTINATION

The TRAVELER PLUS is intended for use indoors or outdoors in pedestrian areas by a single person with reduced mobility up to a maximum weight of 120 kg. Handling is such that it cannot be used by people with severe physical disabilities or balance problems. Your assessor/health professional will determine whether the scooter is suitable for your intended use.

**Note: People with lower body amputations should not use this device unless they are fitted with artificial limbs, as stability may be compromised due to an uneven center of gravity.**

## TECHNICAL SPECIFICATIONS

Class 1 medical device

Overall dimensions	Length 113 x width 57 x height 94 cm
220W motor 2 x 12V 20 Ah batteries Range 20km. Maximum speed 6.6 km/h	
Weight	62 kg
Maximum load	120 kg

\*The weight of the passenger, difficult terrain conditions, low temperatures and the state of the battery may affect the range.

## FEATURES GUIDE

1. Handlebar control panel
2. Removable swivelling padded seats with fold-down backrest
3. Detachable battery pack
4. Non-marking tyres
5. Handlebar angle adjustment
6. LED front lighting
7. Detachable shopping basket

**TO DO :**

- Read this manual carefully before using your scooter.
- Switch off the power supply when transferring between the scooter and your vehicle.
- Choose a slow speed when driving down slopes (max. 6°) or on uneven terrain.
- Keep your feet on the scooter while driving.
- Make sure the seat is in the fixed forward-facing position before driving.
- Make sure the handlebar is securely fastened.
- Make sure the batteries are fully charged before you travel.
- As far as possible, avoid rough or soft terrain and long grass.
- Watch out for other road users and pedestrians.
- Switch off the scooter when not in use.
- Follow the maintenance guide to ensure safe use of the scooter.
- Always remove the keys if you leave the scooter unattended. Use of the scooter by unauthorised persons or children may result in injury.

**DO NOT:**

- Carry passengers.
- Slopes greater than 6°.
- Cross slopes.
- Turn sharply at full speed.

- Get on or off the scooter with the power off.
- Cut the power supply while driving - the scooter will stop abruptly (see emergency procedures).
- Use the scooter where you can walk safely.
- Attempt to walk up or down kerbs.
- Use the scooter if you are under the influence of alcohol or any other medication that may affect lucidity.
- Use a mobile phone while driving.
- Try to negotiate steps or escalators.
- Remain on the scooter while it is being transported in a vehicle.
- Attempt to tow any other device.
- Exceed the stated weight limit.
- Connect any other device to the batteries, as this could cause permanent damage and invalidate your warranty.
- Use batteries/chargers other than those supplied by the manufacturer.
- Travelling on water or soft ground (mud, sand, etc.)
- Lifting the scooter by the bonnet or any other plastic part.
- Shifting your centre of gravity significantly, as this could cause the scooter to become unstable.
- Freewheel on hills. This will damage the motor and brake and invalidate the warranty,
- Tow the scooter with any other vehicle in the event of a breakdown.

Everything has been done to ensure that your scooter provides you with a safe and reliable service, provided that you follow the advice given above.

## ASSEMBLY AND ADJUSTMENT

### ASSEMBLY

Before using your scooter, you need to install the handlebar and seatpost.

To fit the handlebar, you will need the handgrip and spring. Fit the handgrip through the right-hand hinge.

Position the spring on the screw. Align the two hinge joints and place the handle thread in the right-hand side of the hinge. Adjust the handle to the desired angle and tighten the handle fully.

There are two electrical connections - white and black plugs. These must be connected to the corresponding coloured connector. Once connected, carefully push the connectors one by one through the hole in the chassis cover and insert the eyelet.

To fit the seat post, the rear of the scooter must be removed to gain access to the fixing bolt. Locate the seat post and align the holes in the frame and upright at the desired height. Tighten the bolt using the Allen key and the spanner supplied. Refit the rear of the scooter.

For transport, isolate the battery by removing the fuse. It must be replaced before use. The scooter will not operate until this fuse is fitted. The fuse is packed in the small polythene bag with the tools and the seat post bolt.

The fuse must be installed in the fuse holder located behind the charging socket cover at the front of the battery pack.

## **ADJUSTMENTS**

### **Adjusting the angle of the handlebar :**

At the bottom right of the scooter is the angle adjustment knob. To adjust the angle, loosen the knob. Once you have reached the desired angle, tighten it. Make sure the bar is securely fastened before using the scooter.

### **Swivel seat :**

Seat swivel / Seat depth

The seat has a lockable swivel base to make transfers easier. Simply push the lever (circled below) forward and turn the seat at the same time. Releasing the lever locks the seat into position. The seat locks every 45 degrees. Always ensure that the seat is in the locked position before mounting or dismounting the scooter. The seat must always be locked facing forward when driving.

The seat depth can be adjusted by lifting the yellow lever (arrow at bottom) on the lower left-hand side of the seat. Slide the seat into the desired position and release the lever to lock it in place.

### Removing the seat :

To remove the seat, fold the seat down and away from the chassis by pushing the seat swivel lever forward. Do not hang luggage or other objects from the back of the seat, as this could affect its stability.

### Adjusting the seat height :

Remove the seat and the rear of the scooter. Using the tools provided, unscrew the locking bolt from the frame and adjust the seat to the required height. Replace the locking bolt.

Reassemble the rear of the scooter and install the seat.

***The higher the seat is adjusted, the higher the centre of gravity and the more stability may be affected.***

### Adjusting the armrests :

On each side of the seat towards the rear are two knobs and locking pins, which are used to adjust the width of the armrests or remove them for transport.

Turn the knobs anti-clockwise to loosen and remove the locking pins.

This will allow you to adjust the armrests to the desired width. Re-tighten the knobs when you have reached a comfortable width and re-insert the locking pins completely. Remove the armrests completely for transport.

When removing the armrests for transport, make sure the handles are tight to prevent them from coming loose.

Do not apply your full body weight to the armrests when transferring. Lift your arm away from the armrest when transferring to another seat.

## Handlebar control panel :

### A - Battery gauge

This gauge indicates the approximate level of charge remaining in the scooter's batteries. The measurement will be more accurate when the scooter is in motion. **Always fully charge the batteries before driving your scooter and never allow the batteries to discharge completely as this could permanently damage the batteries.**

### B - Horn

The button activates the scooter's horn. **This should only be used to warn others of your presence and not as a reprimand of any kind.**

### C - Maximum speed button

The speed control allows you to set the maximum speed of the scooter. Turn clockwise to increase the maximum speed and anticlockwise to decrease it. **Always start at a low speed until you get used to the characteristics of your scooter. When driving in pedestrian areas such as shopping areas and busy pavements, always respect other pedestrians and reduce the maximum speed so that the scooter moves at the walking speed of others around you.**

### D - Start/Stop button

The key switch is located on the top of the body. Simply insert the key and turn it clockwise to switch on the device. The display lights up to indicate that power is on. The electronic controller checks the scooter's circuits, so there may be a slight delay before the levers are operational. **Remove the ignition key when leaving the scooter parked to prevent unauthorised use. Do not turn off the power to stop the scooter unless it is an emergency. The scooter will stop very suddenly if stopped in this way.**


### E - Front lighting

This button is used to switch on the LED headlight.

## DISASSEMBLY AND REASSEMBLY

The scooter can be quickly and easily dismantled for transport. Always remove the seat and battery before attempting to dismantle the scooter.

Lift the front of the scooter by the seat post while pulling up on the release lever indicated by the arrow below. The scooter will now be in two parts.

To reassemble, allow the rear section of the scooter to rest upside down on the anti-tip wheels. Lift the front section by the seat post towards the rear and place the 'ears' on either side of the frame over the main tube of the rear section as shown below. The front section can now be gently lowered into position.

Once lowered, press down on the seat post and you will hear the front and rear section click into place. To ensure a secure connection, pull on the seat post to ensure that the rear section is correctly connected.

Rear frame main tube

Front frame lugs

## USE

### Clutch release / Freewheel operation:

The clutch release lever is located on the right rear of the scooter. This lever allows you to disengage the drive mechanism and push the scooter into "freewheel" mode. To disengage the drive, push the lever forwards until you reach the stop position (indicated by direction 'N').

Simply push the lever towards the rear of the scooter to re-engage the drive (position 'D'.) The scooter will not operate when the lever is in the disengaged position. The scooter must be switched off and on using the key switch to reset the circuits before normal operation can be resumed. The scooter will beep continuously if the brake is disengaged when it is switched on.

**Caution! Never leave the scooter unattended in freewheel mode, as it may roll down a slope. Freewheel mode should only be used in an emergency to move the scooter due to a breakdown. NEVER move the scooter in this mode when it is occupied. NEVER tow the scooter.**

In freewheel mode, the scooter can be pushed on level ground at walking speed. The scooter will brake automatically if this speed is exceeded.

**NEVER sit on the scooter and coast on a slope as this can cause permanent damage to the motor, brake and controller due to excessive regenerative braking.**

### Getting started

#### Starting to drive :

**Note: You must fully charge the batteries before first use.**

Choose an area with plenty of space.

Mount the scooter by lifting the armrest and sliding onto the seat.

**Make sure the seat is securely locked in the forward-facing position before driving.** Adjust the handlebars and seat if necessary (see 'Adjustments'). **DO NOT use the handlebar to provide full support during transfers.**

The speed and direction of the scooter are controlled by the accelerator lever. The arrows in the image above indicate the direction in which the scooter will move when the lever is operated as shown.

#### **At the start**

Insert the key, switch on the power and first set the cruise control to a low level. The further you pull the lever, the faster the scooter will travel, but only up to the maximum speed that has been selected. Release the lever and the scooter will stop smoothly.

**Warning! DO NOT turn off the power to stop the scooter, as this will cause the transmission to lock and the scooter to stop abruptly, which could result in serious injury. This procedure should only be used in an emergency.**

#### **Steering**

The scooter's steering is light and easy. When starting off, make sure there are no obstacles in the way of the rear wheels. When driving in tight spaces, turn the handlebars in the direction you want to go and apply power gently. The scooter can turn a very tight corner with complete stability at low speeds.

**Caution! NEVER turn sharply at high speed as this can make the scooter unstable. The maximum speed should be set at 50%.**

### Reversing

Before reversing, make sure that the area behind the scooter is clear of any obstacles or hazards. Set the speed to low before reversing. Keep clear of all corners, angles and obstacles.

**Caution! Do not reverse on a main road.**

### Ramps / Slopes

The stability of your scooter depends on factors such as the height of the seat, the angle of the slope, and the height and weight of the user. When approaching a slope, do so directly and not at an angle, and avoid turning whenever possible.

### Driving surfaces

Your scooter will perform well on many surfaces. However, avoid long grass, loose gravel, sand, mud and ice.

### Transfer

Raise the armrest on the side you wish to transfer to. Use a balance aid if necessary. **DO NOT use the handlebar or armrests to support the full weight during transfers.**

### Batteries and battery chargers:

#### Batteries

**Fully charge the batteries (8-12 hours) before using for the first time.**

The TRAVELER PLUS uses two AMP 20Amp, 12V batteries.

The batteries are located in battery pack compartments under the front of the seat. They are easily removed by simply lifting them out using the handle(s) on the top, removing the right-hand side battery pack first). Always make sure that the batteries are correctly repositioned before using them.

Low temperatures, uneven terrain and the user's weight can affect battery performance. Your battery gauge only serves to determine the charge level remaining in your batteries and will give you its best indication when the scooter is in motion.

**Note: Only batteries recommended by the manufacturer MUST be used if they are to be replaced, ONLY an authorised TRAVELER PLUS dealer can carry out battery replacement.**

#### Battery charger

The battery charger is a 6 amp outboard switch mode charger. The battery charging port is located on the front of the battery box (see picture below). Before using your scooter for the first time, make sure you fully charge the batteries, this can take up to 12 hours.

To charge the batteries, plug the charger's 3-pin male plug into the charging sockets - see below. Plug the 3-pin mains plug into a wall socket and switch it on. A steady red LED lights up when the battery charger is connected and charging is in progress. The green LED lights up when the charging cycle is complete. Always switch off the charger before disconnecting the plug from the battery casing.


**Caution! Do not charge batteries for more than 12 hours. Always charge the batteries in a well-ventilated indoor area. You must continue to charge the batteries even if the scooter is not to be used for an extended period (2 to 3 weeks). Failure to do so may result in permanent damage and render the batteries unusable.**

As far as possible, do not interrupt the charging cycle until the batteries are fully charged. Continued interruption could permanently damage the batteries. If the charging cable is damaged, stop using the device immediately and ask your retailer to replace it.

**ONLY USE THE CHARGER SUPPLIED WITH YOUR SCOOTER.** Do not leave the charger connected to the charger port once charging is complete. **ALWAYS** remove the keys from the scooter when charging.

Battery life is affected by the frequency of incomplete recharges. Whenever possible, recharge batteries to their maximum capacity.

### **Sunshine**

On sunny days, it is possible for the surfaces of your scooter to reach high temperatures (like those of motor vehicles). Please take precautions if you leave the scooter in direct sunlight for any length of time. Shade or cover the seat and hand controls to minimize any increase in temperature.

### **Health considerations**

The suitability of this product has been assessed on the basis of the user's current state of health. It is essential that the user regularly monitors their ability to use the scooter safely.

### **Product stability**

Scooters are designed to allow the user a certain degree of flexibility of movement in restricted areas. Driving at top speed and attempting tight turns can create a situation where the product becomes unstable. It is therefore essential that the user pays attention to the advice given in this manual.

For maximum stability, the speed setting should be adjusted up to 50% of the maximum speed available (halfway between the tortoise and the hare or by simply entering the red sector of the speed indicator, depending on the model).

The main factors affecting stability are:

- Maximum speed setting
- Seat height
- User movement
- User manoeuvres
- Slopes or uneven terrain

Please keep your speed at a level where you can stop safely without having to make erratic turns.

Do not load bags, etc. on the back of the seat, as this will affect the centre of gravity and could cause instability.

Do not lean out of the seat, especially when moving it. Avoid turning on slopes.

## Driving in pedestrian areas

As mentioned previously, when using the scooter in pedestrian areas, always lower your maximum speed to that of the pedestrians around you. Remember the stopping distance of your scooter and always leave enough space to stop safely to avoid collisions with pedestrians or other objects.

## CARE AND MAINTENANCE

### General cleaning

The shrouds and seat of the scooter should be cleaned with a damp cloth and a mild detergent. Do not use abrasive cleaners or solvents, as this will damage the plastic components. Do not spray water (from a hose or high-pressure cleaner) onto the scooter, as this could damage the electronics.

### Routine maintenance

We recommend that you have your scooter serviced at least once a year by an authorised dealer. Please check the following:

Check the product	Anytime	Weekly	Monthly
Check freewheel operation	X		
Checking your scooter's speed	X		
Scooter cleanliness			X
Check the condition of the battery			X
Check the condition of the tyre cable sheaths			X
Check that the seat is in the correct position for use	X		
Check the brakes	X		
Check the condition of the tyres		X	

## TROUBLESHOOTING

Your scooter is equipped with the latest electronic controls programmed to protect the electrical system against abnormal overloads.

### Fuses

There are two fuses on the battery cables to protect against possible overloading of the wiring.

### Checklist

If your scooter does not work, please check the following:

- The unit is switched on.
- All connectors and connections are secure.
- The battery charge level indicates that the battery is fully charged.
- The freewheel device is in the driving position.
- The circuit breaker has not tripped if fitted.

If the freewheel device has been used with the scooter switched on, the scooter will not operate until the key switch has been turned off and on again.

If a fault occurs, the fault will be identified by a sequence of beeps which relate to the fault in the table below.

Bip Code	Description	Scooter status	Action
1	Charge the batteries	Scooter can drive slowly	Batteries less than 23.3v. Charge as soon as possible.
2	Poor motor connection	The scooter does not work	Check all connections between the motor and the controller.
3	Motor short circuit	The scooter does not work	Check all connections between motor and battery.
4	Not in use		
5	Not in use		
6	Inhibition circuit activated	The scooter does not work	Check that the charger is not yet connected or that there is a controller fault.
7	Accelerator pot error	The scooter does not work	Check that the wig lever is in the neutral position. Faulty potentiometer.
8	Controller fault	The scooter does not work	Check that all connections are secure. Switch the scooter on and off, this may clear the code.
9	Magnetic brake error	The scooter does not work	Check that the scooter is not freewheeling. Check the brake connections.

**NOTE: If your scooter does not work after checking these possible problems, consult your dealer immediately.**

**If the scooter is switched on and not used for a period of time (approximately 10 minutes), it may enter a "sleep" mode. Simply switch the scooter off and on again to reset the controller.**

#### **Safety information on electromagnetic interference (EMI)**

Scooters and electric wheelchairs are designed to operate under certain conditions. However, radio waves or electromagnetic fields can affect the operation of the scooter.

Radio or television transmitters or portable devices such as mobile phones, portable CB radios, etc. can be a source of electromagnetic interference.

#### **Emergency procedures**

If your scooter behaves erratically or becomes unresponsive while driving, release the handlebar lever so that it returns to the neutral position and turn off the power using the key or switch. This will stop the scooter and activate the braking system, preventing any further uncontrolled manoeuvres.

**NOTE:** Turning off the power using the key or switch may cause the scooter to stop abruptly and should only be used to stop the scooter in an emergency, as this could damage the engine.

## CHANGE OF OWNERSHIP

The product can be reused. It must first be cleaned and disinfected in accordance with hygiene instructions. When transferring the product, remember to hand over all the necessary technical documents to the new user. Beforehand, the product must be inspected by an authorised specialist.

## RECYCLING

When the product becomes unusable and you need to dispose of it: please contact your specialist retailer. If you wish to recycle it yourself, ask your local waste management company for disposal instructions.

### Disposing of and recycling used batteries

**Batteries must always be recycled correctly. Consult your local regulations.**

## STORAGE

**Failure to observe the storage conditions may lead to deterioration of the product and thus to the risk of serious injury.**

Do not store your product for prolonged periods near a source of heat or in direct sunlight (e.g. behind a window or near a radiator) or near a source of cold.

Keep away from open flames and sources of sparks.

Observe the following conditions when storing your rollator:

- In a dry, temperate place
- Protected by packaging, from dust and corrosion (e.g. abrasive elements, sand, seawater, salty air).
- Store all removed parts together, in the same place (or mark them if necessary) to avoid mixing them with other products when reassembling.
- All components must be stored without any load (do not place heavy objects on product parts, do not jam anything...).
- If the scooter gets wet, wipe it with a dry cloth to ensure that no water penetrates the internal components of the scooter.
- For everyday use, do not use in rain showers or puddles.

If you are not going to be using the scooter for a while, we recommend the following:

- Fully charge the battery
- Remove the battery pack
- Store the scooter in a dry place away from extreme temperatures.
- If the battery is stored with a low charge, if you do not charge the battery often, or if you store the scooter in extreme temperatures, the battery may be seriously damaged.

**Warning! Do not charge the battery if it has been stored below 0°C. This could cause serious damage to the battery.**

When storing the scooter for a long period, it is advisable to protect the floor with a blanket or cardboard. This will reduce tyre wear.

## PERIOD OF USE

Your scooter has been designed to give the user a long and reliable life, provided it is properly maintained and regularly serviced by an authorised dealer. However, due to inevitable wear and tear and technological improvements, it is considered that the average usable life of this scooter is five years.

## WARRANTY

This product is guaranteed for two years from the date of purchase by the user. IDENTITIES will provide replacement parts for any part that has broken due to a manufacturing, design or material defect. Contact the distributor from whom you purchased the product in case you need to make a claim. Failures resulting from accidental damage, misuse, alteration, or wear and tear are not covered by this warranty. Contact the dealer from whom you purchased the product to apply the warranty.

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