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INTRODUCTION

Thank you for purchasing this blood pressure monitor. Please read this manual carefully before use. Contact your physician if you have any questions regarding your blood pressure.

This fully automatic device uses the non-invasive oscillometric method, which detects blood movement in your brachial artery, to measure your blood pressure and heart rate. The results are displayed on a digital screen. Without using a stethoscope, you can quickly and easily obtain the measured values.

This device complies with Regulation (EU) 2017/745 on medical devices, as evidenced by the CE marking accompanied by the notified body identification number.

This device complies with the following standards :

- EN ISO 81060 concerning non-invasive sphygmomanometers – Part 1: Requirements and test methods for non-automated measurement devices.
- EN 1060-3 concerning non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems.
- EN 60601-1-2 on medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests.
- EN 1060-4:2004 on non-invasive sphygmomanometers – Part 4: Test procedures to determine the overall accuracy of automated systems.
- ISO 81060-2:2018 concerning clinical validation of automated non-invasive sphygmomanometers.
- IEC 80601-2-30:2018 on medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

Intended use

This device is intended for use by healthcare professionals or by individuals at home. It is designed to measure systolic and diastolic blood pressure in adults using a non-invasive technique, in which an inflatable cuff is wrapped around the arm.

SAFETY WARNINGS

CAUTION indicates a hazardous situation which, if not avoided, could result in serious injury or death.

WARNING indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

CAUTION

- Consult your physician before starting to use this blood pressure monitor.
- As with all oscillometric blood pressure measuring devices, certain conditions may affect measurement accuracy, including:
 - Cardiac arrhythmias
 - Severe hypotension
 - Low blood perfusion
 - Shock
 - Diabetes
 - Vascular abnormalities
 - Individuals with implanted electrical devices (e.g., pacemakers)
 - Pregnant women
 - Preeclampsia
 - Patient movement during measurement
- In such cases, the oscillometric method may produce incorrect values. This poses a health risk if the values are misinterpreted. Always consult your physician to interpret the results.
- This product does not provide a medical diagnosis. The results are for reference only. Self-diagnosis and self-medication based on the measured values may be harmful to your health. Always consult a physician for any treatment or medication adjustments.
- Technical measurement tolerances are possible. See the section “TECHNICAL SPECIFICATIONS”.

- Frequent repetition of measurements may cause side effects, including :
 - Nerve compression with temporary paralysis of the hand
 - Release of an arterial or venous blood clot, which may be life-threatening. Consult your physician for specific risks related to cuff use in your case
- The “irregular heartbeat” function does not replace a medical examination but may help detect certain abnormalities at an early stage. Consult your physician for confirmation.
- This function is not intended to diagnose or treat arrhythmia. Only a physician can do so.
- The WHO blood pressure classification chart is for reference only. It does not replace a medical diagnosis.
- In case of unusual or suspicious variations in measurements, consult a physician immediately.
- Women who have undergone a mastectomy or axillary lymph node removal should consult their physician before using this device.
- The device must be used under the specified environmental conditions, otherwise measurements may be inaccurate. See “TECHNICAL SPECIFICATIONS.”
- Do not apply the cuff anywhere other than the left arm. Incorrect use may be harmful to your health.
- Do not use on children, infants, or individuals unable to give consent (e.g., mental disorders). Consult a professional for alternative methods suitable for children.
- Packaging materials pose a suffocation risk to children. Dispose of them immediately and keep them out of reach.
- This product contains small parts that may be swallowed. Keep out of reach of children.
- Cuff size is crucial. Use only on adults whose arm circumference matches the specified range (see “TECHNICAL SPECIFICATIONS”).
- Electromagnetic interference: Keep the device away from sources of disturbance (e.g., mobile phones, microwaves). This may affect measurement accuracy.
- Batteries must not be recharged or reactivated by other means, as this may cause an explosion.
- In case of battery leakage, avoid any contact with fire, as this may cause ignition or explosion.
- If battery fluid comes into contact with eyes or skin, do not rub. Rinse thoroughly with clean water and consult a physician.
- Do not use the device in the presence of flammable gases (e.g., anesthetic gases, oxygen, hydrogen) or flammable liquids (e.g., alcohol).
- Use only the cuffs and accessories recommended by the manufacturer. Use of non-compatible components may damage the product or endanger your health.
- The tube may pose a strangulation hazard. Keep the device out of reach of children and individuals requiring supervision.
- Do not wrap the tube around the neck. Risk of strangulation.
- Do not use during transport or movement (e.g., car, ambulance, helicopter, walking, exercise), as this may distort the measurement.
- Use, transport, and store the device under the conditions specified in this manual. Otherwise, the results may be inaccurate.


WARNING

- Do not use the device if your arm is injured or if it contains a catheter. This could worsen the injury.
- Remove all jewelry or accessories from the arm before measurement. These could cause bruising.
- Do not place the cuff over thick clothing (e.g., jacket, sweater), as this may distort the measurement and cause hematomas or skin marks.
- If the cuff continues to inflate, immediately stop the measurement by pressing the START/STOP button, then remove the cuff.
- In case of battery leakage, avoid direct contact. Wear gloves and clean the compartment with a dry cloth.
- Do not disassemble the device. This could be dangerous. If a problem occurs, refer to the “TROUBLESHOOTING” guide or contact an authorized service center.
- Ensure there are no folds when applying the cuff to avoid bruising.
- Temporary marks may appear on the skin after a measurement, especially if repeated, in hypertensive patients or those with a weak pulse. In rare cases, these marks may persist for several days. Consult your physician to evaluate the risks in your case.
- Do not place anything on the tube during measurement (e.g., arm or objects), as this may distort the results.

MAINTENANCE AND STORAGE

- The blood pressure monitor consists of high-precision electronic components. The accuracy of the measurements and the service life of the device depend on careful handling. Protect the device from shocks (e.g., falls), moisture, water, dirt, dust, chemicals, extreme temperatures, strong temperature fluctuations, direct sunlight, and nearby heat sources (e.g., stoves, radiators). These may damage the device. It must be stored under the specified environmental conditions. Please refer to the “Technical Specifications” section.
- The device is designed and manufactured for long service life. However, it is generally recommended to have the device inspected every 2 years to ensure proper functioning and measurement accuracy. Contact your dealer for servicing.
- Never immerse the device or its components in water, and do not spill any liquid on it, as this may cause damage.
- Never use rechargeable batteries. This may damage the device.
- Replace all batteries at the same time with batteries of the same type. Do not mix new and used batteries.
- Never attempt to repair, open, disassemble, or adjust the device yourself. This could damage the device or impair its functions. In case of malfunction, follow the instructions in the “Troubleshooting” section or contact your dealer.
- Do not allow any objects to fall into the openings or tubes of the device. This may cause damage.
- Do not press the buttons with excessive force or with sharp objects.
- Clean the device and the cuff with a soft, slightly damp cloth, then immediately dry with a soft, dry cloth. Do not rub.
- Never use harsh solvents, detergents, or strong chemicals (e.g., thinners, alcohol, benzene).
- When storing, do not place any heavy objects on the device.
- Do not excessively bend the cuff or the tube. Avoid sharp folds and keep them away from sharp edges.
- Leaking batteries may damage the device. If you do not plan to use it for an extended period, remove the batteries from the compartment before storing it.
- Used equipment, components, and batteries must not be disposed of with household waste. They must be disposed of in accordance with local regulations. Improper disposal may cause environmental pollution.

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
Button operation [1] / [2]

- Normal measurement with reliable data
→ Press the [1] / [2] button
- Data recall
→ Press and hold the [1] / [2] button
- Unrecorded measurement (guest mode)
→ Press both buttons simultaneously [1] / [2]


MEASUREMENT DISPLAY

- Air release – Libération de l'air
- Pressurization – Pressurisation
- Date – Date
- Time – Heure
- Severe hypertension – Hypertension sévère
- Moderate hypertension – Hypertension modérée
- Mild hypertension – Hypertension légère
- High normal – Haute Normal
- Normal – Normale
- Optimal – Optimal
- Average of previous measurements – Moyenne des mesures précédentes
- Memory – Mémoire
- Systolic pressure – Pression systolique
- Diastolic pressure – Pression diastolique
- Pulse rate – Fréquence du pouls
- Irregular heartbeat indicator – Indicateur de rythme cardiaque irrégulier
- Battery indicator – Indicateur de batterie

INSTALLATION / REPLACEMENT OF BATTERIES

1. Push the cover located under the device as shown.
2. Remove the used batteries and insert new ones. Use LR6 or AA alkaline batteries and only use batteries of the same type.
3. Check the polarity (+ and -) according to the markings in the compartment.
4. Replace the battery cover.


NOTES :

- If the batteries are not inserted correctly, the device may not work or may be damaged.
- The icon [] indicates that the battery power is sufficient.
- When the batteries are low, the icon [] and the code “E6” appear. Replace all batteries. Do not mix new and used batteries.
- Never use rechargeable batteries. This could damage the device.
- Battery life varies depending on ambient temperature. It is shorter at low temperatures.

USING AN OPTIONAL AC ADAPTER (sold separately)

1. Insert the AC adapter cord into the jack located on the right side of the device.
2. Plug the AC adapter into a wall outlet.
3. To remove the AC adapter, first unplug it from the wall outlet, then disconnect the cord from the device.

Note :

- The optional AC adapter must comply with IEC 60601-1 requirements.
- Use only the dedicated AC adapter specified by authorized dealers. Other AC adapters may differ in output voltage and polarity, which could pose a risk to your safety and damage the device.
- When the adapter is connected, the device does not use batteries.

SETTING THE DATE AND TIME

When installing new batteries

1. “YEAR” flashes on the screen.
2. Press the button [] / [] to set the current year.
3. Press the button [] to confirm; “MONTH” will then flash.

4. Press the button [⏸1] / [⏸2] to set the current month.
5. Press the button [⌚] to confirm; “DAY” will then flash.
6. Press the button [⏸1] / [⏸2] to set the current day.
7. Press the button [⌚] to confirm; “HOUR” will then flash.
8. Press the button [⏸1] / [⏸2] to set the current time.
9. Press the button [⌚] to confirm; “MINUTE” will then flash.
10. Press the button [⏸1] / [⏸2] to set the current minutes.
11. Press the button [⌚] to confirm. The date and time setting is complete.

At Any Time

1. Press the button [⌚] to turn on the device when it is off.
2. The date and time will be displayed.
3. Press and hold the button [⌚] for 3 seconds, and “YEAR” will flash.
4. Follow the same procedure as above to set the clock.

TO APPLY THE CUFF

1. Pass the end of the cuff farthest from the tube through the metal ring to form a loop. The smooth fabric should be on the inside of the loop.
2. Insert your left arm into the cuff loop. The bottom edge of the cuff should be positioned about 2–3 cm (1 inch) above the elbow. The white artery marker must be placed over the brachial artery on the inside of the arm. The tube should run down the arm and be aligned with the middle finger. Do not place the cuff over thick clothing (e.g., jacket or sweater), as this will prevent the monitor from working properly. Blood circulation must not be restricted by tight clothing or other objects.
3. Pull the cuff so that the top and bottom edges fit snugly around your arm.


4. When the cuff is positioned correctly, press the sewn hook material firmly against the rough surface of the cuff.
5. Make sure the cuff fits snugly on your arm and is in direct contact with the skin.
6. If the cuff is properly applied, the hook material will be on the outside of the loop, and the metal ring will not touch the skin.
7. Insert the cuff tube into the jack located on the left side of the device. Make sure it is securely inserted and that no bends obstruct the tube.

BEFORE MEASURING BLOOD PRESSURE

1. Apply the cuff according to the instructions in the section *"TO APPLY THE CUFF."*
2. Rest for at least 5 minutes before each measurement; otherwise, the results may be inaccurate.
3. Sit in a comfortable position: keep your legs uncrossed, feet flat on the floor, and your back supported by a chair. Your arm should rest on a flat surface with the cuff at heart level.
4. Relax your arm with the palm facing upward.
5. Stay relaxed, keep still, and do not talk during the measurement to avoid inaccurate results.

NOTES :


- This device is intended for use in adults only. Do not use it on or by children, toddlers, or infants.
- For reliable monitoring, it is recommended to measure your blood pressure at the same time every day.
- To obtain resting blood pressure, do not eat, drink alcohol or caffeinated beverages, smoke, exercise, or take a bath for at least 30 minutes before measurement.
- Blood pressure varies depending on your recent activities. To minimize fluctuations due to physical activity, rest for 5 to 10 minutes before taking a measurement.
- Avoid taking measures if you are stressed or physically fatigued.
- Perform measurements in a quiet, relaxed place at room temperature.
- Do not talk, remain still, and stay relaxed during the measurement.
- Always wait at least 5 minutes between two measurements to allow blood circulation to return to normal.
- If an abnormality is detected, the device will stop the measurement and display an error code. See *"TROUBLESHOOTING."*
- If you have a very weak or irregular heartbeat, the device may have difficulty measuring your blood pressure accurately. (An irregular rhythm is defined as a variation of 25% compared to the average of beats during the measurement.)
- Do not place any object on the tube during the measurement (arm, objects, etc.), as this will distort the results.

TO MEASURE YOUR BLOOD PRESSURE

Measurement with data storage

1. Follow the steps in the section "BEFORE MEASURING BLOOD PRESSURE."
2. Press any of the buttons [] / [] when the device is off to start the measurement. The current pressure will be displayed as the measurement begins.
3. The cuff will start to inflate. It is normal for it to tighten firmly around the arm. A pressure bar will appear during the measurement. See "PRESSURE INDICATOR" for more details.
4. Once inflation is complete, the cuff will automatically deflate. When the pulse is detected, [] will flash with each beat: the measurement is in progress.
5. At the end of the measurement, the systolic pressure, diastolic pressure, and pulse rate will be displayed and stored. The cuff will then fully deflate.

NOTE :


You can stop inflation at any time by pressing any of the following buttons [] / []

Measurement without Data Storage (Guest Mode)

1. Follow the steps in the section "BEFORE MEASURING BLOOD PRESSURE."
2. Press any of the buttons [] / [] while the device is off to start the measurement. The average blood pressure and the last data number will briefly appear. The current pressure is displayed when the measurement begins.
3. The cuff will start inflating. It is normal for it to tighten firmly around the arm. A pressure bar will be displayed during the measurement. See "PRESSURE INDICATOR" for details.
4. Once inflation is complete, the cuff will automatically deflate. When the pulse is detected, [] will flash with each beat: the measurement is in progress.
5. At the end of the measurement, the systolic pressure, diastolic pressure, and pulse rate are displayed and recorded. The cuff then fully deflates.


NOTE :

You can stop inflation at any time by pressing any of the buttons [] / []. The measurement result will not be stored in the average data or memory.

PRESSURE INDICATOR


This indicator monitors the progress of the pressure during measurement.

WHO CLASSIFICATION INDICATOR

Each segment of the bar indicator corresponds to a blood pressure classification according to the WHO.

WHO Classification :


TO CHECK THE AVERAGE AND PREVIOUS MEASUREMENT DATA

This device is equipped with a memory function that records the measurement results for each user. At the end of each measurement, the result is automatically stored.

1. To obtain the average data of the last 3 measurements and display the previous measurement data, press and hold one of the user buttons [81] / [82] for 5 seconds. The average data will be displayed.
2. Press the same user button again and hold it continuously to display the previous measurement data.
3. Press and hold the same user button for 5 seconds to turn off the device.

TO DELETE MEASUREMENT DATA

1. Simultaneously press the buttons [81] / [82] and hold them down when the device displays the average or previous measurement data of the selected user.
2. [CL] and [--] will appear. Keep holding the buttons until [CL] and [00] are displayed.
3. All measurement data of the selected user will be deleted.

WHAT IS AN IRREGULAR HEARTBEAT ?

This blood pressure monitor can measure blood pressure and pulse rate even in the presence of an irregular heartbeat. A heartbeat is considered irregular when it varies by 25% from the average of the beats recorded during measurement. It is essential to remain relaxed, still, and silent during the measurement.

NOTES :

- We recommend consulting your physician if the [♥] indicator appears frequently.
- The “Irregular Heartbeat” function is not a substitute for a medical examination but may help to detect potential irregularities at an early stage. Always consult your physician to determine what is appropriate for your case.
- This function is not designed to diagnose or treat arrhythmia. Only a licensed physician can diagnose arrhythmia.

ABOUT BLOOD PRESSURE

What is blood pressure ?

Blood pressure is the force exerted by the blood against the walls of the arteries. Systolic pressure corresponds to the contraction of the heart. Diastolic pressure corresponds to its relaxation phase. Blood pressure is expressed in millimeters of mercury (mmHg). An individual's natural blood pressure is represented by the baseline pressure, measured upon waking, at rest, and before any food intake.

What is hypertension and how can it be controlled ?

Hypertension is abnormally high blood pressure which, if left untreated, can lead to various health problems, including strokes or heart attacks. It can be controlled through lifestyle modifications, stress reduction, and medical treatment under medical supervision. To prevent or manage hypertension :


- Do not smoke
- Engage in regular physical activity
- Reduce salt and fat intake
- Have regular medical checkups
- Maintain a healthy body weight

Why measure blood pressure at home ?

Measuring blood pressure in a medical office can cause stress and distort the results (an increase of 25 to 30 mmHg compared to a home measurement). Home measurement reduces the influence of external factors, complements the readings taken by your physician, and provides a more complete and representative history of your blood pressure.

WHO Blood Pressure Classification

Evaluation standards for hypertension, regardless of age, have been established by the World Health Organization (WHO).

TROUBLESHOOTING

Nothing appears on the screen, even when the device is switched on	The batteries are exhausted	Replace all batteries with new ones.
	The battery polarities are not positioned correctly	Reinstall the batteries, ensuring that the (+/-) polarities match the indications in the compartment.
	The plug or contact at the power outlet is loose (if using an AC adapter)	Check the wiring to ensure that the plug and power supply are properly connected.
Error Code 1 (E1)	The cuff is incorrectly positioned	Sit comfortably, remain still, and ensure the cuff is at heart level.
Error Code 2 (E2)	You moved your arm or body during the measurement	Remain perfectly still and silent during the measurement.
Error Code 3 (E3)	The cuff is not attached correctly. The cuff cannot be applied.	Attach the cuff properly. Check that the cuff tube connection is securely attached to the device.
Error Code 4 (E4)	The device is unable to measure	If your pulse is too weak or irregular, the device may have difficulty determining your blood pressure.
	Measurement error	Sit comfortably, remain still, and reposition the cuff correctly.
Error Code 5 (E5)	Excessive cuff inflation	The measurement exceeds the 300 mmHg limit. It is recommended to consult a physician as soon as possible.
Error Code 6 (E6)	Low battery level	Remaining power is insufficient for operation. Replace the batteries with new ones.
The blood pressure monitor continues to inflate	Circuit blockage	Remove and reinstall the batteries, then repeat the measurement.

TECHNICAL SPECIFICATIONS

Measuring Range	Systolic Pressure: 50 – 250 mmHg Diastolic Pressure: 30 – 200 mmHg Pulse: 40 – 180 beats/minute
Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of the displayed value
Resolution	Pressure: 1 mmHg Pulse: 1 beat/minute
Measuring Method	Non-invasive oscillometric method
Power Supply	4 $\times$ 1.5 V AA batteries (approx. 500 uses)
AC Adapter (optional)	Input: 100 – 240 V Frequency: 50 – 60 Hz Output: 6 V @ 600 mA (service life: approx. 5000 uses)

Operating Temperature / Humidity	+5 °C to +40 °C 15 to 90% RH (non-condensing)
Storage Temperature / Humidity	–25 °C to +70 °C Up to 90% RH (non-condensing)
Atmospheric Pressure for Operation, Storage, and Transport	700 hPa to 1060 hPa
External Dimensions	Approx. 99 × 160 × 56 mm
Arm Circumference	22 – 36 cm (Standard), 17 – 22 cm (Optional), 32 – 44 cm (Optional)
Accessories	Cuff, User Manual, Storage Pouch, Batteries (optional), AC Adapter (optional)
Classification	Applied part: Type BF  Equipment: Class II 

Remarks :

- Mode of Operation: Continuous
- Caution: Refer to the SUPPLEMENTARY DOCUMENTS

APPENDIX I

Recommendations and Manufacturer’s Declaration – Electromagnetic Emissions			
The blood pressure monitor (MD2520/MD2540/MD2580) is intended for use in the electromagnetic environment specified below. The customer or user of the blood pressure monitor (MD2520/MD2540/MD2580) should ensure that it is used in such an environment.			
Conducted and Radiated RF Emissions	Compliance		Electromagnetic Environment – Recommendations
RF Emissions CISPR 11	Group 1		The blood pressure monitor (MD2520/MD2540/MD2580) uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause interference with nearby electronic equipment.
RF Emissions CISPR 11	Class B		The blood pressure monitor (MD2520/MD2540/MD2580) is suitable for use in all facilities, including domestic environments and those directly connected to the public low-voltage power supply network intended for residential use.
Recommendations and Manufacturer’s Declaration – Electromagnetic Immunity			
The blood pressure monitor (MD2520/MD2540/MD2580) is designed to be used in the electromagnetic environment specified below. The customer or user of the blood pressure monitor (MD2520/MD2540/MD2580) should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Recommendations
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power Frequency (50 Hz – 60 Hz) Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.

Recommendations and Manufacturer's Declaration – Electromagnetic Immunity

The blood pressure monitor (MD2520/MD2540/MD2580) is intended for use in the electromagnetic environment specified below. The customer or user of the blood pressure monitor (MD2520/MD2540/MD2580) should ensure that the device is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V/ms 150 kHz to 80 MHz	3 V	Mobile and portable RF communication equipment should not be used in close proximity to any part of the blood pressure monitor (MD2520/MD2540/MD2580), including the cables, at a separation distance recommended and calculated from the equation applicable to the transmitter frequency.
Radiated RF electromagnetic fields IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power of the transmitter (in Watts), d is the minimum recommended separation distance (in meters), and E is the immunity test level (in V/m). Electromagnetic fields emitted by fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level for each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol : 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

The field strengths of fixed transmitters, such as base stations for radio telephones (cellular/cordless) and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted with accuracy in theory. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the blood pressure monitor (MD2520/MD2540/MD2580) is used exceeds the applicable RF compliance level above, the blood pressure monitor (MD2520/MD2540/MD2580) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the blood pressure monitor (MD2520/MD2540/MD2580).

Within the frequency range of 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

APPENDIX III

Recommended separation distances between portable and mobile RF communication equipment and the blood pressure monitor (MD2520/MD2540/MD2580)

The blood pressure monitor (MD2520/MD2540/MD2580) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the blood pressure monitor (MD2520/MD2540/MD2580) can help prevent electromagnetic interference by maintaining a minimum distance between portable or mobile RF communication equipment (transmitters) and the blood pressure monitor (MD2520/MD2540/MD2580), as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter (W)	Separation distance according to transmitter frequency (m)		
	150 kHz to 80MHz D=3.5VP	80 kHz to 800MHz D=1.2VP	800 kHz to 2.5GHz D=2.3VP
0.01	0.117	0.12	0.23
0.1	0.370	0.38	0.73
1	1.170	1.2	2.3
10	3.700	3.8	7.3
100	11.7	12	23

For transmitters rated for maximum output power not listed in the table above, the recommended separation distance d (in meters) can be estimated using the equation applicable to the transmitter frequency, where P is the transmitter's maximum output power in watts (W), as specified by the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the applicable separation distance is that of the higher frequency range.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

APPENDIX IV

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception (which can be determined by turning the equipment off and on), the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

WARNING

- To comply with the limits for Class B digital devices pursuant to Part 15 of the FCC Rules, this device complies with Class B limits. All peripherals must be shielded and grounded. Use of uncertified peripherals or unshielded cables may cause interference to radio or television reception.
- Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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Symbol for the marking of electrical and electronics devices according to Directive 2012/19/UE.
The device, accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal.

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
Grandway Technology (Shenzhen) Limited
No. 5, the Second Industrial Zone, Zhukeng Community, Longtian Street, Pingshan District, Shenzhen,
518118 The People's Republic of China

Shanghai International Holding Corp. GmbH (Europe) Eiffestrasse 80, 20537 Hamburg, Germany.

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