COUSSINS GRAINES DE LIN NUQUE – UNIVERSEL – LOMBAIRE – LONG LOMBAIRE/CERVICAL - CHAUSSONS

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EN

PRECAUTIONS FOR USE

PLEASE READ CAREFULLY THE PRECAUTIONS FOR USE AND INSTRUCTIONS BEFORE EACH USE!

- Be careful when removing the cushion from the microwave, as it may be very hot.
- To avoid frostbite or burns, let the cushion sit until it reaches a skin-tolerable temperature.
- Apply for no more than 20 minutes at a time, with 20-minute breaks between applications.
- To prevent fire hazards, DO NOT exceed the maximum heating time.
- DO NOT apply hot or frozen beads directly to your skin (use the plush side).
- DO NOT wear while sleeping.
- Check your skin frequently during use.
- DO NOT use on unhealthy, damaged, or broken skin, or on bruised or swollen areas that have occurred within the last 48 hours, or on areas where you do not feel heat or other sensations properly.
- DO NOT use on individuals unable to remove the product themselves, including children, infants, the elderly, or anyone unable to follow all usage instructions.
- Do not overheat.
- DO NOT apply additional pressure to the product.
- DO NOT ingest.
- Keep out of reach of children and pets.
- DO NOT use if the heat pack's contents are leaking and/or the cover is damaged or torn.
- If the contents come into contact with your eyes, remove the cushion immediately, rinse the affected area with water, and seek medical attention.
- Discontinue use and consult a physician if pain worsens or persists for more than 7 days, or if you experience discomfort, burning, swelling, rash, or any changes in your skin.
- Consult a doctor before use if you have diabetes, poor circulation, heart disease, rheumatoid arthritis, or if you are pregnant.

USE

1. Wrap the cushion in a moist towel before heating. To prevent burns, make sure the product is fully wrapped in a

damp towel or cloth. Flip the fabric each time before reheating. Place the cushion at room temperature, evenly spread out in the center of the clean microwave turntable, without touching the microwave walls.

2. Make sure all additional microwave settings, such as "browning" mode, are turned off.
3. Heat in the microwave for the time specified in the table below.
4. Check that the temperature is appropriate before applying. If additional heat is needed, reheat for the time indicated below. Before reheating, check if any part of the cushion is too hot. Knead and rotate the cushion to distribute the heat evenly. Do not overheat.

ITEM	MAXIMUM CONTINUOUS HEATING TIME (DO NOT EXCEED)	INITIAL HEATING FROM ROOM TEMPERATURE	REHEATING TIME IF MORE HEAT IS NEEDED
835193 – Neck Cushion 835195 – Lumbar Belt 835197 - Slippers	90 seconds (800W) 80 seconds (1000W)	60 seconds (800W) 50 seconds (1000W)	10 seconds (800W) 10 seconds (1000W)
835194 – Universal Cushion	110 seconds (800W) 90 seconds (1000W)	80 seconds (800W) 60 seconds (1000W)	10 seconds (800W) 10 seconds (1000W)
835196 – Long Lumbar/Cervical Cushion	150 seconds (800W) 120 seconds (1000W)	120 seconds (800W) 90 seconds (1000W)	10 seconds (800W) 10 seconds (1000W)

CARE INSTRUCTIONS

Spot clean only. Lay flat to dry. Do not bleach. Do not machine wash, iron, or tumble dry. Do not immerse in water. May be stored in the freezer in a sealed bag until use, or kept at room temperature.

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