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INSTRUCTIONS MANUAL EN

IMPORTANT

Read the instructions carefully before using this product and keep them in case of need.

User: If you are unable to read or understand the warnings, precautions or instructions, please contact your healthcare professional or distributor before using this product, as this may cause risks and damage. They will be able to advise and inform you on the most appropriate use. Consult a physician before using any medical device. He or she will advise and inform you on the most appropriate use of this product

Distributor: this leaflet must be given to users of this product.

PRODUCT DESTINATION

Memory-effect repositioning aid with re-ventilation, intended for bedridden patients

This product is intended for bedridden and immobile patients, for the prevention of pressure ulcers. It works by reducing pressure points on the head. This device is intended for healthcare professionals and requires a medical prescription. Read the properties of the cushion carefully before use.

PRECAUTIONS FOR USE

This product must be used on the advice of a healthcare professional following a medical prescription.

If the patient does not accept the product, do not use it. If the patient does not tolerate the product, do not use it.

If the patient experiences pain when using this product, do not use it.

Follow the instructions in the medical prescription and the planned treatment protocol.

Do not use a soiled or undecontaminated cushion.

Do not use the product if it is damaged - Check the seams before use.

Failure to observe the precautions for use may result in discomfort and injury to the hemiplegic patient.

If the patient does not accept the product, do not use it.

USE

It is advisable to position the patient every 2 hours in a 30° supine position (semi-fowler) and in a 30° left and right lateral position. Read each cushion specification and possible combinations. Ask a healthcare professional for advice.

Clinical monitoring involves ensuring that the patient is properly hydrated and eating properly.

CHANGE OF PATIENT

When the product is to be re-used on another patient: it must be decontaminated and used in perfectly hygienic conditions.

When changing patient, remember to give this leaflet to the new user.

CARE AND MAINTENANCE

The STERITEX cover for the POSITPRO cushion can be decontaminated with a non-chlorinated disinfectant (spray, wipe). The cover should be machine-washed at 95°C with the zip closed. Disinfect the cushion after each use. Dry with a soft, clean, dry cloth. Do not use abrasive products, detergents or chlorinated products. In the event of wear or tear, do not attempt to repair the product. Dispose of it or contact your retailer to find out if it is covered by the after-sales service.

STERITEX cover

Cushion core

RECYCLING

When the product has become unusable and you need to dispose of it, contact a local waste disposal company.

If you wish to recycle the product, please enquire locally.

STORAGE

Failure to comply with storage conditions may lead to deterioration of the product and therefore a risk of injury.

The product can be stored between 0°C and 50°C.

Keep away from sources of ignition.

Observe the following conditions when storing your products

- In a dry, temperate place
- Protect your product with dust-proof packaging
- Do not place heavy objects on the cushion: the foam could settle and deteriorate.

TECHNICAL CHARACTERISTICS

Assembly: polyester-cotton thread stitching

Filling: Comfortable TECHNI GEL fibre, exceptionally soft and supple to the touch. 100% polyester with memory effect + foam pad made of viscoelastic.

Core cover : Made from 100% cotton warp and weft fabric.

Cover in contact with the patient: Made from STERITEX technical textile, polyurethane coating, soft and comfortable to the touch (190g/m²).

STERITEX polyurethane is breathable and waterproof thanks to its 100% polyester bi-elastic knitted base. (190g/m² - Composition: 55% polyester - 45% polyurethane)

This product contains biocides in accordance with EU Standard 528/2012. Protection against microbiological deterioration.

Active ingredients: silverglass complex 398477-47-9 + zinc - pyrethione

The covers are fitted with a ZIP with flap for better protection of the cushion.

Overall dimensions	829154 – 45 x 40 x 8 cm
Weight	829154 – 0.900 kg

PERIOD OF USE

The expected life of this product is 2 years under normal use and maintenance. After this period, if the product has lost its memory properties, if the cover is damaged or if the Velcro fasteners are damaged: do not use this product. It can be used as long as it is still in good condition.

The lifespan of the product will vary according to the weight of the patients who use it and the number of times it is washed. The flattening of the product will be an indicator of its age.

WARRANTY

This product is guaranteed for two years from the date of purchase by the user against defects in stitching and workmanship. Contact the distributor from whom you purchased the product should you need to make a claim. Failures resulting from accidental damage, misuse, improper storage, improper unpacking, modifications or wear and tear are not covered by this warranty.

Contact the retailer from whom you purchased the product to apply the guarantee.

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